



Laser Safety Advice and Equipment for Laser & IPL Clinics

If you are running a health, beauty, cosmetic or tattoo clinic and are using lasers or IPL for skin resurfacing, hair removal, tattoo removal, physiotherapy, dental work or similar treatments; Lasermet are the right people to help you with laser safety compliance and with the Healthcare Commission laser regulations.

We specialise in laser & IPL safety equipment and advice. Our aim is to help you deal with all your laser & IPL safety requirements as quickly and easily as possible, thus allowing you to get on with running your business in the knowledge that your staff and patients are fully protected from any laser hazards. In order to enable this we supply all the laser safety related products and services which you are likely to need, including:

- Healthcare Commission Registration
- Laser Protection Adviser Service
- Laser Safety Training for your Staff and Clinicians
- Laser Protective Eyewear for Laser Operator
- Laser Protective Eyewear for Patients
- Laser Blocking Curtains and Screens
- Illuminated Signs
- Laser Interlock Systems
- Labelling and Other Signage
- Laser Fume Extraction

Healthcare Commission Registration

If all the paperwork involved in registering your clinic with the Healthcare Commission for laser and IPL treatments is causing you a headache, let us organise it for you. Our experienced LPA's will sort out your registration quickly and efficiently, allowing you to get on with running your clinic.

Laser Protective Eyewear for Laser Operator

We supply a very large range of laser protective eyewear with a great choice of frame styles. See our eyewear brochures, browse the website or call our sales staff for advice on appropriate eyewear.



Laser Protection Adviser Service

All clinics using Lasers or IPL systems are required by the Healthcare Commission to have a Laser Protection Adviser (LPA). Lasermet LPA's are accepted by the Healthcare Commission under the Association of Laser Safety Professionals' Certification Scheme. Where required we will construct your local rules, provide a complete Laser Safety Management File and make sure that you are complying with HC requirements.

IPL Protective Eyewear

The intense light from IPL systems can cause eye damage. We offer a range of eyewear for operators and also for patients (see below).



Laser Safety Training for your Staff and Clinicians

The Healthcare Commission require that all clinicians and staff who are present during operation of your laser or IPL system must be proficient in the Core of Knowledge defined in the Guidance on the Safe use of Lasers in Medical & Dental Practice. We supply in house training, tailored to your specific needs and led by lecturers highly experienced in the subject.

**67 Portchester Road
Bournemouth
BH8 8JX United Kingdom
Tel: +44 (0) 1202 770740
Fax: +44 (0) 1202 770730
e-mail: sales@lasermet.com
Website: www.lasermet.com**

Laser Protective Eyewear for Patients

A range of eyewear for patients designed to provide maximum protection and combined with maximum access to the face.



Laser Blocking Curtains and Screens

Lasers can be dangerous at large distances from the source. Laser blocking roller blinds or curtains are often necessary to stop a dangerous laser beam passing through window. Room divider systems are also available.



Labels, Illuminated Signs & Laser Interlock Systems

We offer a range of illuminated signs which will warn people not to enter the room during operation of the laser. For maximum protection for staff and patients alike use our laser interlock control system with automatic switching of illuminated signs. We also supply labelling and other useful signage.



Laser Fume Extraction

Fume extraction is necessary to protect the operator from viral infections which may be carried in the plume. Filtration is also desirable to remove bad odours produced by vapourised hair etc.



Lasermet offer an extensive range of laser safety products and services.

LASERMET LTD
laser safety solutions 